**MANUAL DE USUARIO  
USER'S MANUAL  
MANUEL DE L'UTILISATEUR  
BETRIEBSANLEITUNG**

**Cabinas estériles por flujo laminar horizontal  
Horizontal flow clean benches  
Cabine stérile de flux laminaire horizontal  
Steril-kabinen mit laminarer horizontalströmung**

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CABINAS HORIZONTALES  
MODELOS: AH, MH, BH

**USER'S MANUAL**

*HORIZONTAL FLOW BENCHES  
MODELS: AH, MH, BH*

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## **GB** 1.- GENERAL DATA

*Operating Instructions Manual for horizontal laminar flow sterile benches, which ensures high level protection for the sample:*

Models: - **AH-100**  
- **MH-100**  
- **BH-100**

### **WARNING**

*These operating instructions must be read before starting up the laminar flow bench.*

*Follow all the safety measures given in these operating instructions.*

*Unplug the bench before attempting any Maintenance operations.*

*Any use or application of the bench not according to the recommendations stated in these operating instructions may originate a dangerous situation and is the cause of the loss of the manufacturing guarantee.*

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## **GB** 2. BENCH DATA

### 2.1 Operation

The horizontal flow benches are equipped with high efficiency centrifugal fans, assuring high protection for the sample.

The outer safety cabinet is made of oven dried lacquered sheet steel, functional in design, small outside dimensions and front access for filter changing.

The air is blown by the fans (1) into chamber/plenum (2) where it is filtered through the HEPA absolute filter (3) and then arrives at the working area (4) in laminar form; this process obtains Class 10 according to the US. Fed. Std. 209 D (see Annex A1.).

### 2.2 Application

This model of horizontal laminar flow sterile bench is designed for handling samples or products when the Class 10 sterile protection must be guaranteed in the working area, provided that no protection to the user or the environment is required.

See principal diagram in Annex A1.

### **WARNING**

Never handle products or samples in presence of corrosive, inflammable or explosive gasses.

Never handle products or samples biologically pathogenic or likely to harm in any way the user or the environment.

Never handle cytostatics (chemotherapy) or radiological samples.

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See Paragraph 10.

**2.4 Certification**

Control of quality and certification with the results of the trials carried out according to Norms. The certificate is located in the interior of the safety cabinet.

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### **GB 3. TRANSPORT**

The bench is shipped in single enclosure and in such a way to avoid transport damage.

**This Manual as well as the control certificate are supplied with the safety cabinet.**

#### **4. INSTALLING INSTRUCTIONS**

For a correct horizontal laminar flow bench start-up procedure, follow these instructions step by step.

##### **4.1 Placement**

- a) Cabinets should be located away from doorways, passageways, air diffusers and the like, which would influence cabinet airflows.
- b) The safety cabinet should be located firmly on the table or bench so that no vibrations are produced.
- c) It is advisable to keep the room temperature within the limits of working comfort in order to prevent over-heating in the working area.
- d) Modification of room air inlet diffusers to deflect air away from cabinets or a reduction in room inlet air velocity should be considered where necessary.  
  
Sidewall air diffusers are unsuitable for rooms where cabinets are installed.
- e) Personnel traffic generates air movements, which could cause a sudden loss of barrier containment. Warning signs or other measures to inhibit traffic during critical manipulations should be considered.

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**GB** *Partitioning of cabinet locations may be used to reduce the effects of air movement from other cabinets, and from personnel traffic.*

- f) *All windows should be fixed.*
- g) *The sudden opening and closing of doors can disturb cabinet airflows, and should be avoided when cabinets are in use.*
- h) *The cabinet location will directly influence the sound level. The level will increase significantly in small rooms with highly reflective surfaces, because of the proximity of cabinets to walls and where fan speeds have been increased to compensate for room air currents or increased filter resistance.*
- i) *Cabinets should be re-tested prior to use. Changes made in the room or the air supply system could significantly change airflow characteristics.*

#### 4.2 Warnings for the first connection

- *Before connecting to the mains, clean the bench in a general way so that any dust particles accumulated during transport are removed. Follow the cleaning and disinfecting instructions.*
  - *Make the required electrical connection (see the description plate).*
  - *Before working for the first time, start the bench to blow through the filters and leave it on for about 6 hours.*
- The new filter gives off a characteristic smell, which tends to disappear after a short while.*

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### **GB 4.3 Electrical connection**

The safety cabinet will be supplied with single-phase electric fans, connection cable (on request), general switch and two protecting fuses.

The flow switch is located on the control panel (Fig. 2 pos. 1).

The electronic control panel has protecting fuses.

#### **WARNING**

Before switching on the bench check that the supply mains coincides with the data given in the characteristics plate.

The bench may be activated only with the switch placed in the top right side.

The client's plant or machinery where the bench is connected must comply with the prevailing safety rules in every country.

#### **4.4 Starting up**

After dealing with the above points, turn the key switch (located at top-right of the bench) to start up the laminar flow bench. Ensure that the suction air intake is clear.

To switch the fans on, just push "FAN" (Fig 2. pos.1) in the control panel.

#### **4.5 Switching off**

To stop the bench, just push "FAN" (Fig. 2 pos. 1) in the control panel.

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## **GB 4.6 Stopping for long periods**

Turn on the general switch on the top of the cabinet, just close to the power cable.

In the event of a long stop, avoid as much as possible the accumulation of dust. When using the bench again, proceed as if it were started for the first time. See paragraph 4.2.-

## **5. INSTRUCTIONS FOR USE**

### **5.1 General recommendations**

- The bench working area must not be used as laboratory equipment storage. This may cause an unnecessary presence of dust, which endangers sterility.
- Even if the working area is in a clean sterile laminar flow, the bench surroundings may be contaminated. The contaminating sources must be detected in order to avoid them while introducing arms, material, etc.
- All the material necessary for the work must be particle free and previously cleaned.
- Do not put inside the working area materials such as: paper, wood, cardboard, pencils, erasers, etc. since they give off a lot of particles.
- Is recommended to wash arms, hands and finger nails with a germicidal soap before and after working. The users must avoid touching their mouth and eyes.
- Is recommended to use long sleeved overalls with tight cuffs, and, with special work, protecting gloves. (Alternatively use mittens). Both the overall and the mittens must be made of a material, which gives off as few fibres and particles as possible.

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- GB** - Start the bench 10 minutes before beginning your work with it. Thus the working area and the material introduced are scavenged and the particles removed.
- Whenever pipettes are used, these must be of mechanical suction. Never use them with mouth suction since you are very likely to inhale the aerosol originated by the suction.
  - Whenever platinum handles must be used, we suggest the use of electrical incinerators and even disposable ones, if possible.
  - If the work to be done requires the use of a Bunsen burner, or similar gas flame, we recommend those with a constant pressure button. Remember that an excessively long flame might burn the absolute filters and bend the light diffuser.
  - The shadows and whirls caused by the objects, equipment and material surrounding the working area must be examined before starting the operation, in order to determine their possible effects. Remember that the airflow does not go back to its laminar condition before a distance 2.5 times the diameter of the object causing the obstruction.
  - When vials and tubes are used, we recommend the use of screw caps instead of cotton ones, the latter gives off a lot of particles.
  - Once the operation is finished, all disposable products (handles, Petri Plates, etc.), and the culture medium, samples, tubes, vials, etc. will be removed from the bench inside impermeable and sterilizable bags, if necessary.
  - When working with the bench or during its cleaning operations, special care must be taken not to damage the HEPA filters by hitting them, pouring liquids, splashing, etc.

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
**GB 5.1.1 Ergonomics**

For their own security, the user should keep in mind the positions adopted during the job carried out in the safety cabinet.

Find as follows, examples of positions NOT recommended

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
Besides as protection so much the product as of the user, is necessary the safety cabinet to be located in a place where any current can be produced, that can cause turbulence in the laminar flow.

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## GB 5.2. Control panel

### Components description (Fig. 2).

2x16 LCD display with 3x5.5 mm characters for the different indications and programs.

Six buttons.

Control dial functions:

- Button (1): operates the fan (on/off). By default, the fan speed is set as standard, if you wish it to go slower, you should then press the down button (nr 6).
- Button (2): operates the light (on/off), linked electrically with the UV germicidal lamp button (nr 3).
- Button (3): operates the UV germicidal lamp (on/off) (optional), linked electrically with the lighting switch (nr 2).
- To connect the UV lamp with a timer, you should press the button (nr 3) and start it up. The button (5) activates or deactivate the timer.
- Steps to follow to program the timer for the UV lamp:
  - 1.- You should start with the UV lamp switched off.
  - 2.- Press the down button (nº6).
  - 3.- Select the time, from 1 to 59 minutes, using the up and down buttons (nr 5/nr 6).
  - 4.- Confirm with the UV's start button (nr 3).
- Button (4): operates the AUX output for possible safety valves, service plugs, etc.

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### GB Signal lights, alarms and display:

- Signal light (nr.7) for air-impulsion system functioning, green.
- Signal light (nr.8) for lighting switched on, blue.
- Signal light (nr.9) for UV germicidal lamp switched on, yellow.
- OPTIONAL. Filter clogging alarm (nr.11), red.
- UV germicidal lamp on (nr.12), red.
- Display functions:


Fig.2

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- Kit lámpara germicida (U.V.) (\*)

- GB**
- 1- The display usually shows TELSTAR and the version of the program.
  - 2- Viewer for hours of air impulsion system and UV germicidal lamp functioning V. (VENT/FAN: 000000h, U V: 000000h) will be activated with button 5, only when the rest of the functions are stopped.
  - 3- Normal or slow fan speed (velocity) (FAN VEL. STD or ALARM FAN VEL 1/2).
  - 4- Slow speed air impulsion alarm system (ALARM VEL. 1/2).
  - 5- Fix or programmed period of time for UV germicidal lamp (ALARM UV ON or UV ON 12 MIN).
  - 6- UV lamp alarm on (ALARM U. V. ON).
  - 7- Filter clogging alarm (ALARM FILTER).

### 5.3 Optional accessories

The horizontal laminar flow benches can have several additional accessories according to the work requirements.

These optional accessories are the following:

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## GB WARNING

*When this accessory is used, please take the usual precautions with this type of lamp. It should only be used for the previous external sterilisation of the elements introduced into the bench. Therefore, this lamp is an accessory previous to the operation with the bench, and it is essential to have it turned off when the user is in front of the table or working at it. The SECURIT glass panes protect from direct radiation but the UV rays are reflected on bright metal parts and might damage the worker's eyes. This is why the germicidal lamp must only be used when nobody is working with the bench.*

- Gas stop-cock/tap (\*)
- Vacuum stop-cock/tap (\*)
- Other gasses stop-cock/tap (nitrogen, etc.) (\*)

## WARNING

*When compressed air or nitrogen are used, use an absolute filter in the line to avoid the introduction of particles inside the working area.*

- Support table (\*)
- Six hook fluid therapy rod (\*)
- HEPA/ULPA absolute filters clogging optical alarm.

**(\*) The user referring to the diagrams added to these Operating Instructions can easily assemble these accessories. Fig. 4, 5, 6 and 7.**

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## GB 6. MAINTENANCE INSTRUCTIONS

### WARNING

Unplug the bench before attempting any Maintenance or disassembly work. Usually one regular test should be performed every year to check all the elements of the bench.

### 6.1 Maintenance Table

The frequency given in the maintenance table for each operation, and especially that of the filter replacement, are recommended indicative values for normal use of the bench. Harder working conditions may shorten the intervals stated in the following table.

Maintenance table

OPERATION	PERIODICITY
Filter leakage test	Once a year minimum
Downflow velocity	Once a year minimum
Smoke visualisation of product protection	Once a year minimum
Filters pressure drop test	Once a year minimum
Replacing absolute impulsion and exhaust filters	Every 4000/5000 working hours, depending on quality of ambient air

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## GB 6.2 Regular validation and control

The International Regulations recommend regular validation and control (at least once a year). The tests to be performed must be the same as those done at the last manufacturing stage, issuing a result of certification.

### 6.3 Prefilters

The horizontal laminar flow sterile benches have some prefilters in the air intake section (Fig. 8 pos. 3). These must be cleaned or vacuum cleaned every 500 hours and replaced every 1000 hours.

#### WARNING

Never start up the bench without having previously installed the prefilter correctly, since this would cause an early HEPA filter clogging.

### 6.4 Absolute filter

The average life of the HEPA filters with these benches is approximately 4000/5000 hours, which are logged and indicated by the machine's elapsed time meter (Fig. 2 display).

After this period of time, the blown air volume decreases due to the higher filter resistance, and therefore results in lower bench performance.

#### WARNING

In order to be able to guarantee the right assembling and fitting of the filter, as well as the lack of leaks, the filter replacement operations should be performed by specialised technicians, with the necessary checking tools: air speed indicator, leak detector, aerosol generator, etc.

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
**GB 7. BENCH ASSEMBLY AND DISASSEMBLY INSTRUCTIONS**

- 1- Fit the head piece (2) and fasten it to the frame (1) with the bolts M8.
- 2- Fit the inlet screen (3) over the head piece (2), tighten with the bolts M5.
- 3- Fit the 2 glass side walls (4).
- 4- Disassembly is the reverse of the assembly procedure.

**Fig.3**

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## GB 8. INSTRUCTIONS TO ASSEMBLE OPTIONAL ACCESSORIES

### 8.1 U.V. Germicidal kit

The Germicidal Kit is easy to install, following the instructions (Fig. 4). The "UV" button on the control panel (Fig. 2 pos. 3) switches it on. This button is electrically linked with the light button for safety reasons.


#### WARNING

The UV lamp radiation is harmful to people. Absolutely never work while the U.V. lamps are on.

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#### 8.1.1 Instructions

- Remove fluorescent lamp from germicidal kit.
- Connect plugs 6 and 7.
- Put the germicidal kit in the bottom of the fan front panel and screw the bolts M5 in the 4a-4b and 5a-5b holes.
- Put the fluorescent lamp back into its place.

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
Fig.5

## GB 8.2 Stop cock/Tap for gas, vacuum, air, nitrogen, etc.

On the right hand side, at the back, there are two holes to house the intakes for these fluids; follow the instructions. (Fig. 5).

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- Remove the glass cap (1) and the adhesive (2) from the bench.
- Place the stop-cock (3) in the side hole.
- Place the rubber joint (4) and the washer (5).
- Fasten the stop-cock (3) using the pipe nipple and the nut (6).

## F 8.2 Cannette/robinet gaz, vide, air, azote, etc.

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**GB**


**8.3 Support table**

All the models of horizontal laminar flow benches may be placed on top of a support table, which is very easy to install (Fig. 6).

**8.3.1 Instructions**

1- Support leg	5- Hexagonal screw M6x10, washer M6
2- Support union legs	6- Sunk screw M6x50, washer M6
3- Support front	7- Rubber support AR-M10, hexagonal nut M10
4- Cover	

- a) Fasten the support union legs (2) on the support legs (1), using the washers and hexagonal screws M6 (5).
- b) Put the support front (3).

**Fig. 6**

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
## GB 8.4 Fluid therapy rod

Stainless steel bar with six hooks to hang vials or bags for fluid therapy (Fig. 7).

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- Disassemble the aluminium blade (1) and make two Ø 10 holes corresponding with the cabinet serblocks M6 (8) and put it back in its place.
- Screw-on the rod supports (5) in the serblocks M6 (8).
- Shift the side wall glass forward slightly (4).
- Introduce the stainless steel rod (2) through the first support, put the hooks (6) and introduce it through the last support.
- Put the side wall glass back to its place (4).


Fig. 7

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
**Fig. 2**

**Fig. 3**

## **GB ANNEX A2. CLEANING AND DISINFECTION**

For worktable and sides, use sterile or disposable cloths, which do not shed particles or fibres.

To disinfect use cloths dampened with a disinfecting solution which will not damage or affect the paintwork, stainless steel or glass. Cleaning prior to disinfecting can also be done with water and soap.

To access corners or hard to reach areas, other types of cleaning tools can be used (as shown in the Figures 2 and 3).

For inside the safety cabinet, cleaning is recommended on the following occasions:

- Before starting any work in the safety cabinet.
- After working in the safety cabinet.
- Whenever there is a change of work programme.
- In event of liquid spilling on the work table.
- Before carrying out mechanical or biological control tests in the working area.

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## GB ANNEX A6. IMPULSION FILTERS INTEGRITY TEST

### PURPOSE

Check the integrity of the downflow filters.

### SCOPE

AH, MH and BH.

### TOOLS

- Aerosol generator (cool or hot).
- Shell Ondina 15 or Emery 3004.
- Photometer (ATI TDA 2-G).

### METHODOLOGY

- Turn on the photometer and adjust it for measurements in accordance with the manufacturer's instructions (Fig. 1).
- Inject aerosol into the flow before the filter (Fig. 2).
- With the nozzle of the probe approximately 2 cm. from the surface, scan the downstream side of the HEPA filters (including the perimeter of each filter pack) by passing the photometer probe in slightly overlapping strokes over the entire surface.

Also scan the entire periphery of the filter and the junction between the filter to filter mounting frame at a scanning rate that does not exceed 5 cm/s (Fig. 3).

Fig. 2

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**GB ACCEPTANCE CRITERIA**


It is considered a "leak" when the photometer reading in the area behind the filter detects a concentration greater than 0.01% of that obtained in the part in front of the filter.

**Fig. 3**

**Figura 4.** Resultado correcto: no fuga  
*Correct result: filter not leaking*  
*Résultat correct: pas de fuit*  
*Korrektes Ergebnis: Undicht*

**Figura 5.** Resultado incorrecto: si fuga  
*Wrong result: filter leaking*  
*Résultat incorrect: fuite*  
*Nicht korrektes Ergebnis: Dicht*

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**GB ANNEX A6.1 DOWNFLOW AIR VELOCITY TEST**
**PURPOSE**

To check that the speed control in different points of the working area is suitable to achieve product protection. To determine the speed average as well as the uniformity range in the laminar flow zone.

**SCOPE**

AH, MH and BH.

**TOOLS**

The tool to use is an anemometer

**METHODOLOGY**

The test must be done where the equipment is working. Depending on the model of the cabinet, take as many readings as indicated in the following table:

**ACCEPTANCE CRITERIA**

The average of the airflow velocity at the work surface should be 0,40 m/s ± 20%.

Modelo / Model Modèle / Modelle	Dimensiones / Dimensions Dimensions / Messungen	Puntos / Points Points / Punkte
AH	1220 x 762 x 68	6
MH	1525 x 762 x 68	6
BH	1830 x 762 x 68	8

La sonda del anemómetro se colocará a una distancia de 200 ± 50 mm de la superficie filtrante, en las cabinas de flujo laminar según figuras 1 y 2.  
The anemometer's probe is placed at a distance of 200±50 mm from the filtering surface in the laminar flow cabinets, as in figures 1 and 2.  
Effectuer les lectures 5 à 10 cm au-dessus de la vitre d'ouverture frontale (selon Fig. 1 et 2).  
Die Messungen sind 5 bis 10 cm oberhalb der Glasscheibe der Frontöffnung durchzuführen.

Fig. 1


Fig. 2

Se calcula el promedio mediante la siguiente fórmula:  $V_{Promedio} (m/s) = 1/n \times (V1 + V2 + V3 + \dots + Vn)$ .  
For this calculation, we proceed to the following formula:  $Average (m/s) = 1/n \times (V1 + V2 + V3 + \dots + Vn)$ .  
La moyenne est calculée grâce à la formule suivante:  $VMoyenne (m/s) = 1/n \times (V1 + V2 + V3 + \dots + Vn)$ .  
Der Durchschnitt wird mit Hilfe der folgenden Formel gebildet:  $Mittlere G. (m/s) = 1/n \times (V1 + V2 + V3 + \dots + Vn)$ .

FILTOS FILTER FILTRÉS FILTER V (m/s)	PUNTOS / POINTS POINTS / PUNKTE							
	1	2	3	4	5	6	7	8
Velocidad promedio / Average Velocity Vitesse moyenne/ Mittlere Geschwindigkeit (m/s)								

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**GB ANNEX A6.3 LIGHT TEST**
**OBJECT**

To check that the level and uniformity of the lighting system is correct, in order to guarantee a safety working area inside the cabinet.

**SCOPE**

AH, MH and BH

**TOOLS**

Light meter.

**METHODOLOGY**

An imaginary division of the working surface will be done in equal parts (not more than 8, and not less than 6).

The illumination level will be measured in each of these parts, by placing the light meter over the working area.

The average of the readings obtained will be calculated, as well as the individual safety illumination levels of the safety cabinet.

**ACCEPTANCE CRITERIA**

The test will be correct if:

- The illumination level average is not less than 750 lux.
- Any individual reading can be less than the 70% of the average illumination level.

Punto / Point/ Point / Punkt	1	2	3	4	5	6	7	8	Promedio / Average/ Moyenne / Durchschnitt
Valor / Value/ Valeur / Wert (lux)									

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Fig. 1

**GB ANNEX A6.5 PRESSURE TEST**
**OBJECT**

To verify the degree of clogging in the filter.

**SCOPE**

AH, MH and BH.

**TOOLS**

Pressure indicator

**METHODOLOGY**

Based on the premise that the initial differential pressure is known (when the filter is new). To measure the differential pressure, proceed as follows:

- Disconnect the pneumatic installation (Fig. 1).
- The pressure indicator is connected (Fig. 2).

**ACCEPTANCE CRITERIA**

If the pressure obtained exceeds the initial differential pressure by between 3 and 4 mmca, the filter is considered to be clogged.

Guidelines for clogging are as follows:

	Differential pressure(mmca)	Differential pressure (Pa)
New filter	10-12	100-120
Clogged filter	14 a 16	140-160

Fig. 2

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## **GUARANTEE CERTIFICATE**

TELSTAR TECHNOLOGIES, S.L. guarantees the substitution of all the defect pieces in the period of TWO YEARS, starting in the date of the expedition of the equipment.

This equipment has undergone satisfactorily all the tests and revisions. Therefore, any damage caused by transport, wrong installation or inappropriate use is excluded of this guarantee.

This guarantee does not include any of the fungible materials comprised in the equipment, such as filters, fluids, filaments, etc.

The inappropriate use of the equipment excludes TELSTAR TECHNOLOGIES, S.L. of any responsibility for eventual direct or indirect damages of any nature.

Likewise, this guarantee will be void if pieces other than the ones supplied by TELSTAR TECHNOLOGIES, S.L. are used for the repairing of the equipment.

**NOTE. This guarantee is valid for all the equipment unless otherwise stated. In order to obtain any information about this guarantee, please contact with our Sales Department.**

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