

# Compact dual path-dual focus SS-OCT system for whole-anterior segment and retinal imaging



Ana Rodríguez-Aramendía<sup>1,2</sup>, Fernando Díaz-Doutón<sup>2</sup>, Jaume Pujol<sup>2</sup>, José Luis Güell<sup>1</sup>, Ireneusz Grulkowski<sup>3</sup>

(1) Instituto de Microcirugía Ocular, Barcelona, Spain.

(2) Center for Sensors, Instruments and Systems Development, Universitat Politècnica de Catalunya, Terrassa, Spain.

(3) Faculty of Physics, Astronomy and Informatics, Nicolaus Copernicus University in Torun, Poland.



## Introduction


Swept source optical coherence tomography (SS-OCT) opens the possibility of imaging the whole eye. Different strategies have been proposed in this line, such as dual path configurations exploiting coherence revival [1] [2], tunable lens technology [3] or interlaced detection [4].

In this work we present a compact dual path-dual focus SS-OCT system that can image sequentially both the whole anterior segment, from corneal apex to posterior face of the lens, and the retina with a FOV of 21°, featuring a novel configuration.

## Methods

### Scheme of the Anterior / Posterior segment SS-OCT

The proposed configuration consists of a dual path-dual focus sample arm that switches between anterior segment and retinal imaging modalities employing a flip mirror. The system is composed by three spectral channels: OCT signal, anterior segment and retinal (1050 nm), eye tracking (850 nm) and vision channel (vis.).




Swept Source: Short cavity swept laser, 50 kHz, 1050nm (Axsun)

Detector: PDB481C-AC (Thorlabs)



Digitizer: 12 bit, 1GS/s (ATS9371, Alazar)

### Optical performance

The optical interfaces have been designed and optimized with optical simulation software (Zemax, Radiant Zemax, LLC)




### Characterization



## Results

A group of volunteers have been imaged with the system in a laboratory environment to assess its performance.



## Conclusions

We have developed a SS-OCT capable of imaging the whole anterior segment and the retina with a dual path configuration. The first measurements on a group of volunteers have shown the potential of the system. A full clinical validation, including eyes with retinal and anterior segment pathologies, will be conducted at IMO.

## References

- [1] McNabb, R.P., Polans, J., Keller, B., Jackson-Atogi, M., James, C.L., Vann, R.R., Izatt, J.A., & Kuo, A.N. *Biomed. Opt. Express* 10, 338-355 (2019).
- [2] Nankivil, Derek et al. *Biomed. Opt. Express* vol. 6,11 4516-28. (2015)
- [3] Grulkowski, I., Manzanera, S., Cwiklinski, L., Sobczuk, F., Karnowski, K., Artal, P. et al., *Optica* 5, 52-59 (2018)
- [4] Kim, Hyung-Jin et al. *Biomed. Opt. Express* vol. 9,2 360-372. (2018)

### Acknowledgements:

ITN BE-OPTICAL H2020-675512  
www.beoptical.eu

